


FULL PACKAGE DOOPANG Items



1 DOOPANG Body



2 Hydrogel Patch




3 Patch Storage Pad



4 USB Charging Cable



Link with **DOOPANG App** for
Smart Migraine Management



A'ali Town, Kingdom of Bahrain
Road 4270 | Block 742 | Bldg. 2504

Phone: +973 17 664521
Mobile: +973 39 606219

E-mail: info@b1group.me
www.b1group.me



“Accessible Mental Health for Everyone”

Medical Device Approved
by the **Korean Ministry of Food and Drug Safety!**

DOOPANG

Medical Device for Relieving Migraine Symptom

Do Not Hesitate to Put It On!

Migraine cannot be completely cured
but could be significantly reduced or controlled.
Thus, early treatment is important.



“Do Not Hesitate to
Put It On”



DOOPANG

Medical Device for Relieving Migraine Symptom - DOOPANG

Manage Your Migraines

Say Goodbye to Migraines via Smart Habits



* Picture of DOOPANG

- ✓ Those who have been absent from school or work due to migraines
- ✓ Those who have difficulties in their academic or work efficiency due to migraines
- ✓ Those who have experienced irritable feelings or melancholy in their daily lives due to migraines



Relieving Migraine Symptom
Reduction in Frequency of Migraine Occurrence



Approval from the Korean Ministry of Food and Drug Safety
Safe Medical Device



Convenient Functionality at a Touch of One Button
Easy Operation




Transdermal Nerve Stimulation (TENS)
Non-drug/Non-invasive

TECHNOLOGY

Learn how DOOPANG works in managing migraine pain, and find out why you need to control your migraine with DOOPANG

The Principle of TENS Technology Introduced by DOOPANG




By applying the principle of transcutaneous electrical nerve stimulation (TENS), DOOPANG delivers electrical stimulation to the upper branch of the trigeminal nerve through the forehead.

By transmitting an electric current from the DOOPANG electrode surface to the skin, action potential is generated in the nerve, thereby relieving migraines and reducing the frequency of migraine occurrence.

58% Improvement in Migraine Pain (when Using Acute Mode)

58% Improvement in Headache Pain after 1 hour and 24 hours



*Source: Acute migraine therapy with external trigeminal neurostimulation (ACME): A randomized, controlled trial/ 2019, Vol. 39, International Headache Society

DUAL MODE

Two Types of Usage Settings, Depending on the Situation

It can be used in dual mode with varying duration and intensity of use.



Acute Mode



It relieves migraine pain.
100Hz
The session runs for 61 minutes.



Prevention Mode



It reduces the frequency of migration occurrence.
60Hz
The session runs for 21 minutes.

EASY TO USE

Anytime, Anywhere with One-Touch Button!

- 1 Wipe the forehead clean.
* If moisture and oil are not completely removed, the adhesive strength of the patch may decrease.
- 2 Apply the patch 1 cm above the eyebrows.
* The patch can be used up to 10 times.
- 3 Set the mode, and start the session in a lying or comfortable position.

